

ARTIFICIAL SEAWATER BASED AMP1 MEDIUM

Table 2. Sea Salt and nutrient additions

<u>Nutrient</u>	<u>Manufacturer/Grade</u>	<u>Primary Stock (M)</u>	<u>Dilution Factor</u>	<u>Final Conc.</u>
NaCl	(2) / SigmaUltra	n/a	n/a	481 mM
CaCl ₂ ·2 H ₂ O	(2) / SigmaUltra	n/a	n/a	10 mM
KCl	(2) / SigmaUltra	n/a	n/a	9 mM
Mg SO ₄ ·7 H ₂ O	(2) / SigmaUltra	n/a	n/a	28 mM
MgCl ₂ ·6 H ₂ O	(2) / SigmaUltra	n/a	n/a	27 mM
NaH ₂ PO ₄	(2) / SigmaUltra	0.025	1:500	50 µM
(NH ₄) ₂ SO ₄	(2) / SigmaUltra	0.4	1:1000	400 µM
NaHCO ₃	(2) / SigmaUltra	0.6	1:100	6 mM
HEPES	(2) / SigmaUltra	0.5	1:500	1 mM
Na ₂ EDTA·2H ₂ O	(2) / 99%	0.012	1:10 ⁵	0.1170 µM
FeCl ₃ ·6 H ₂ O	(1) / Analytic	0.012	1:10 ⁵	0.1180 µM
ZnSO ₄ ·7 H ₂ O	(3) / >99.5%	0.080	1:10 ⁸	0.0008 µM
CoCl ₂ ·6 H ₂ O	(1) / Analytic	0.050	1:10 ⁸	0.0005 µM
MnCl ₂ ·4 H ₂ O	(1) / Analytic	0.900	1:10 ⁸	0.0090 µM
Na ₂ MoO ₄ ·2 H ₂ O	(2) / ACS	0.030	1:10 ⁸	0.0003 µM
Na ₂ SeO ₃	(2) / ~98%	0.100	1:10 ⁸	0.0010 µM
NiCl ₂ ·6 H ₂ O	(1) / Analytic	0.100	1:10 ⁸	0.0010 µM

Manufacturer Index: (1) Mallinckrodt, (2) Sigma, (3) Fluka

NOTE: To avoid contamination, dry nutrients and metals should be weighed without using spatulas.

1. Prepare nutrient stocks as following using Sigma Ultra grade chemicals:

A. 25mM NaH₂PO₄ (pH 7.5)

- Weigh out 0.300g NaH₂PO₄ on dust-free weigh paper
- Transfer into 100mL volumetric flask filled with about 60mL Milli-Q water
- Dissolve NaH₂PO₄ by inverting flask several times
- Adjust to pH 7.5 with NaOH (Mallinckrodt #7708)
- Adjust volume to 100mL mark with Milli-Q water
- Filter with a sterile 0.2µm nylon filtration unit
- Store sterile stock at 4°C

B. 0.4M (NH₄)₂SO₄

- Weigh out 5.284g (NH₄)₂SO₄ on dust-free weigh paper
- Transfer into 100mL volumetric flask filled with about 60mL Milli-Q water
- Dissolve (NH₄)₂SO₄ by inverting flask several times
- Adjust volume to 100mL mark with Milli-Q water
- Filter with a sterile 0.2µm nylon filtration unit
- Store sterile stock at 4°C

C. 0.5M HEPES (pH 7.5)

- Weigh out 59.575g HEPES on dust-free weigh paper
- Transfer into 1L glass beaker with stir bar and 300mL Milli-Q water
- Dissolve HEPES by stirring
- Fill remaining volume to 500mL mark with Milli-Q water
- Check pH by pipetting out small amount onto pH paper for rough estimate
- Turn stirrer on and **SLOWLY** add NaOH, check pH by pipetting out small amount onto pH paper for rough estimate periodically
- When pH is around 7.5 use a calibrated pH meter probe that has been Milli-Q water rinsed and soaked several times (add NaOH as needed to obtain pH 7.5)
- Transfer to volumetric flask or graduated cylinder and adjust volume to 500mL
- Filter, using a 500mL sterile filter system (Corning code# 430773) with .22µm nylon filter and polystyrene bottle (pre-wash the filter 2-3 times with 100mL Milli-Q water to remove residual carbon) attached to a vacuum pump (use 1 inch of Hg)
- Store sterile stock at 4°C

NOTE: We have migrated to using a 3.75mM TAPS buffer as it produces less H₂O₂, but if working with high density cell concentrations then HEPES will work.

D. Working trace metal stock

Use the working trace metal stock from the Natural Seawater Based PRO99 medium recipe, **but ten times LESS of the stock will be added to AMP1 compared to PRO99.**

2. Prepare “Turks Island Salt Mix”³ by adding the following:

dissolve each ingredient sequentially in 2/3 final volume of Milli-Q water

Milli-Q water	1.0L
NaCl	28.10g
Mg SO ₄ ·7 H ₂ O	6.90g
MgCl ₂ ·6 H ₂ O	5.49g
KCl	0.67g
CaCl ₂ ·2 H ₂ O	1.47g

- Fill to appropriate volume with Milli-Q water
- Autoclave and allow to cool

³ The Merck Index (1989) Merck&Co, Inc. (Rahway, N.J.)

3. To Turks Island Salt Mix add the following nutrients stocks:

	<u>1L</u>
NaH ₂ PO ₄ ·H ₂ O	2mL
(NH ₄) ₂ SO ₄	1mL
HEPES	2mL
Working trace metal stock	10µl

4. Store this ASW basal medium at room temperature for up to one month.

5. Prepare 0.6M NaHCO₃ stock fresh on day of medium inoculation (**do not store**)

- Weigh out 2.52g NaHCO₃ on dust-free weigh paper
- Transfer into 50mL volumetric flask filled with 30mL Milli-Q water
- Dissolve NaHCO₃ by inverting flask several times
- Adjust volume to 50mL mark with Milli-Q water
- Filter with a sterile 0.2µm nylon filtration unit

6. NaHCO₃ addition to basal medium:

10 mL of 0.6M NaHCO₃ per 1L basal medium

